



# Rotavapor® R-210

## Superior evaporation

The R-210 represents highest reliability and convenient operation. A wide range of vacuum pumps, recirculating chillers, glassware and accessories enable highest efficiency.

### Reliable

Accurate results  
and high longevity



### Comfortable

Easy operation thanks to the  
electric lift and quick-release  
heating bath



### Customizable

Compatible with a wide range of  
glass assemblies and accessories

## R-210: Your most important benefits



### Reliable

- Accurate results at best performance
- High longevity due to the use of high quality materials (e.g. P+G coated glassware)
- Covers applicative needs with a wide temperature range of the heating bath (20 – 180 °C)
- No thermal overloading in case of power failure due to automatic flask raising



### Comfortable

- Comfortable operation thanks to ergonomic electric lift with quick-action jack
- Easy flask and vapor duct removal thanks to the multifunctional Combi-Clip
- Easy refilling and cleaning due to heating bath with quick-release-system
- Direct access to frequently used settings for intuitive handling



### Customizable

- Precise and reproducible results guaranteed with an optional vacuum controller
- Largest range of applications thanks to seven available glass assemblies
- Different heating baths allow flask sizes up to 5 liters
- Covering most needs due to wide range of complementary products and accessories

## Complete your portfolio



**Vacuum Pump  
V-700**  
The economical  
vacuum source



**Vacuum Controller  
V-850 / V-855**  
Achieve best process  
performance



**Recirculating Chiller  
F-108 / F-114**  
The efficient way of  
cooling



**Multivapor  
P-6 / P-12**  
Efficient evaporation  
for multiple samples

